

Operational Principle of MicroJet Cryo Trap (MJT-1030E)

US Patent: US 6,190,613 B1

Analysis of dilute components in gases or volatile components desorbing from heated sample in a wide band range requires trapping of these components at the head of column by cooling, followed by rapid thermal desorption. We have developed MicroJet Cryo Trap (MJT-1030E) which consists of a mechanism to liquify N₂ gas and micro jet tube for cooling and thermal desorption at the head of column.¹⁾

The flow scheme of MJT-1030E is shown in Fig. 1 and the operation of micro jet tube in cooling and thermal desorption modes is shown in Fig. 2. N₂ gas from N₂ cylinder is liquified through thermal exchange coil immersed in liquid N₂ and is fed to micro jet tube located in the GC oven via N₂ transfer tube. Inside the micro jet tube are a temperature sensor and column holder (metal tube), and the separation column runs through the column holder and is secured. Liquid N₂ jet is blown against the column in the column holder, and the portion of the column is cooled down to -180° C or below. Liquid N₂ goes out from the both ends of the micro jet tube and is rapidly evaporated, preventing moisture to get into the micro jet tube to become frozen. Once the liquid N₂ micro jet is turned off, heated air in the GC oven will rapidly heat the cooled portion of the column at 800° C/min, causing the trapped components to be thermally desorbed.




Fig. 1 Flow Scheme of MJT-1030E

Fig. 2 How MicroJet Cryo Trap Works in Cooling and Desorption modes

1) Hosaka, et al., 3rd Polymer Analysis Symposium, I-6, p15-16 (1998)

Keyword : MicroJet Cryo Trap, Operational Principle

Applications : General Polymer Analysis

Please forward your inquiries via our web page at : <http://www.frontier-lab.com/>, or send us a fax message.

R&D and manufactured by :
Frontier Laboratories Ltd.

1-8-14 Saikon, Koriyama,
Fukushima 963-8862 JAPAN
Phone: (81)24-935-5100 Fax: (81)24-935-5102

Double-Shot Pyrolyzer® is a registered trademark of Frontier Laboratories Ltd.



Scientific Instruments Manufacturer GmbH
Im Erlengrund 21-23
D-46149 Oberhausen
Phone: +49-208-941078-0 Fax: +49-208-941078-88
<http://www.sim-gmbh.de> info@sim-gmbh.de