

<b>Specifications Agilent 1200 Series Preparative Pump</b>	
<b>Flow precision</b>	<0.3 % RSD in range from 1 mL/min to 100 mL/min
<b>Hydraulic system</b>	Dual pistons in parallel
<b>Flow range</b>	Setpoints from 0.001 to 100 mL/min, in 0.001 mL/min increments. Recommended: 5-100 mL/min
<b>Pressure</b>	Operating range 2–40 MPa (20–400 bar, 290-5880 psi) system pressure
<b>Gradient formation</b>	High pressure using two combined pumps
<b>Delay volume</b>	700 µL
<b>pH range</b>	1.0-12.5 recommended